



MIP

The Softer Side of Healthcare



COMFORTABLE SLEEP



SKIN HEALTH



DIGNITY

BARIATRIC SOLUTIONS



WELLNESS



PERSPIRATION CONTROL



CAREGIVER RELIEF

Specialized products designed to meet the needs of your bariatric residents or patients. Our products help provide a healthy skin environment through fabrics with enhanced moisture management properties. Our selection also features thin profile products that are wrinkle resistant to reduce the risk of pressure injuries, commonly experienced by care recipients with limited mobility in bed.

Discover our bariatric solutions under MIP's different brands.

INTEGRIDERM®

SOLUS™
THE PERFECT BALANCE OF COMFORT & DIGNITY

UltraKnit™

FUSION®

Swift



**BARIATRIC SPECIALIZED
PRODUCTS AND ADAPTED
SYSTEMS FOR RESIDENTS
AND PATIENTS**



PERMAPRINT Clothing Protector

A clothing protector specially designed for the larger resident/patient featuring multiple snap positions.

Benefits

- Generous size that offers proper coverage and dignity

Swift **Repositioning System**

The most cost-effective resident/patient repositioning system. Swift UltraSlide helps promote independence and reduces staff/nurse neck and back injuries caused by lifting and repositioning.

Swift Slider (A)

Benefits

- Low-friction backing allows for safe and easy repositioning

UltraSlide Bottom Sheets (B)

Benefits

- Breathable taffeta backing avoids heat buildup
- Additional friction reduction and hence less effort



FUSION **Patient Care Underpad**

By eliminating the potential for pressure points, FUSION creates a smooth surface for a patient/resident's sensitive skin to lay on. With skin integrity at the forefront of its mandate, FUSION helps protect the person's skin from moisture, shear, and friction.

Benefits

- STAY DRY unquilted surface provides superior moisture management





SOLUS™ I.V. Gowns
THE PERFECT BALANCE OF COMFORT & DIGNITY

Solus gowns have been developed using the latest in textile technology to enhance resident/patient care & comfort.

Benefits

- Soft brushed fabric that feels comfortable on the skin
- Generous size that offers proper coverage and dignity
- Moisture management properties keep skin feeling dry

INTEGRIDERM™ Mattresses

Integriderm mattresses offer excellent pressure redistribution for your resident/patient. They are designed to provide a healthy environment for the skin & body.

Benefits

- Breathable and moisture wicking material for better comfort
- 3 Foam layers that provide comfort and support



UltraKnit™ Sheets

Ultraknit sheets are designed with our special 4 way stretch knit fabric that ensures a wrinkle free surface and allows for the total emersion of the resident/patient into the support mattress.

Benefits

- Breathable knit fabric helps control moisture
- Soft texture fabric promotes comfort



MIP BARIATRIC PRODUCT SPECIFICATIONS



PERMAPRINT Clothing Protector

Style	Dimensions	Features
1102443/SKS	24" x 43"	Snap closure

SOLUS™ I.V. Gowns

Styles	Colors/Pattern	Full sweep	Length	Sleeve	Features
PG/BAI-3XL	GREY Nightingale Print	72"	45"	16.5"	Ties and Plastic snaps
PG/BAI-10XL	BEIGE Nightingale Print	107"	46"	19"	

Swift® Repositioning System

Swift® Slider (A)

Style	Backing color	Dimensions
PTD-56/BA/Y	Yellow	50"W x 72"L

UltraSlide Bottom Sheets (B)

Styles	Bindings	Threads	Dimensions
PTD-LS-S-BA-Y	Green	Green	54"W x 84"L x 8"D
PTD-LS-S-BA1-Y	Yellow	Yellow	48"W x 80"L x 7"D

UltraKnit® Sheets

Types	Styles	Bindings	Colors	Dimensions
Bottom Sheet	SK40B/BE	White	Bone	54"W x 80"L x 8"D
Bottom Sheet	SK41B/BE	Orange	Bone	48"W x 80"L x 7"D
Bottom Sheet	SK415B/BE	White	Bone	42"W x 80"L x 7"D
Bottom Sheet	SK42B/BE	Stone Blue	Bone	39"W x 80"L x 7"D
Bottom Sheet * Optimum	SK4200B	Green with yellow treads	White	42"W x 88"L x 8"D
Top Sheet	SK40T/BE	White	Bone	54"W x 80"L x 8"D
Pillow Case	SK1019P/BE	N/A	Bone	19"W x 32"L

INTEGRIDERM® Mattresses

Style	Cover Color	Dimensions
MIP/BAR/SIZE	Navy	On demand*

* Bariatric mattresses are available in 39", 48", 54" or 80" standard mattress widths. Additional sizes are available upon request.

FUSION® Patient Care Underpad

Styles	Colors	Dimensions	Soaker
9118/FP	Blue	34" x 45"	7 oz/y ²
9118 (Only available in US)	Blue	34" x 45"	7 oz/y ²

CONTACT US

mipinc.com

CANADA 1-877-356-2987

USA 1-800-361-4964



WEBSITE



Oeko-Tex Certification

OEKO-TEX.COM