



FUSION UNDERPADS

Patient care underpads

2021 v.1

The Facts

- **Less Deterioration:** Maintains its integrity and performance after industrial washing.
- **No Pilling:** Synthetic surface won't allow pilling to occur.
- **Minimal shrinkage:** Unique laminated construction maintains its size after repeated washes.

Only 3% shrinkage

- **Thin Profile:** Extreme comfort is achieved by the ultra thin profile that lays flat and unobtrusive
- **Eco-friendly alternative to disposable incontinence pads**



FUSION is a square, fused, unquilted underpad that offers the best in reusable bed linen protection.

Specifications

STYLE #	FORMAT	DIMENSION	WEIGHT	SURFACE COLOR	BARRIER
9100	Standard	35" x 35"	16.4 oz	Turquoise	Blue PU
9118	Bariatric	34" x 45"	19.4 oz	Turquoise	Blue PU
9105/FP	Chair pad	17" x 24"	5.4 oz	Turquoise	Blue PU



Contact your sales consultant

mipinc.com



MIP

The Softer Side of Healthcare