



# Microfiber Isolation Gown

2021 v.1

## The Facts:

- Reversible style and stain resistant properties to reduce labour time and costs
- Durable and quick-dry including knit cuffs
- Fabric meets level 2 requirements
- Carbon yarns reduce static build-up



STYLE	FULL SWEEP	SLEEVE INSEAM	WEIGHT	DESIGN
EHX4200-WAT	64"	24"	280 gr	Wrap-Around Coverage with Traceability Grid



Contact your sales consultant

[mipinc.com](http://mipinc.com)



# MIP

The Softer Side of Healthcare