



SYNTHETIC UNDERPAD

2021 v.1

The Facts

- Bonded quilted synthetic underpad with lightly brushed barrier
- 100% polyester, napped white Sonoma surface
- Waved quilted surface provides increased comfort and improved durability
- Needled non-woven soaker
- High performance at value price



Specifications



Contact your sales consultant

mipinc.com



| STYLE # | DIMENSION | WEIGHT | BARRIER | SOAKER |
|----------|-----------|---------|-----------|---------------------|
| 11026/LV | 34" x 36" | 15.3 oz | Dark Blue | 7 oz/y ² |

MIP

The Softer Side of Healthcare